

Bio-Rad Laboratories  
Contact person: Jason Lai  
E-mail address: jason\_lai@bio-rad.com  
Phone: (510) 741-6170  
Address: 4000 Alfred Noble Dr  
Hercules, CA 94547  
Web site: www.bio-rad.com/diagnostics

1. What is the name of your company's point-of-care testing product or product line?  
(For multiple products/lines, please complete NO MORE THAN THREE separate questionnaires.)

TOX/See™ Rapid Urine Drug Screen Tests

2. What does this test detect?

TOX/See™ detects Amphetamine (1000 ng/mL), Barbiturates (300 ng/mL), Benzodiazepines (300 ng/mL), Buprenorphine (10 ng/mL), Cocaine (300 ng/mL), Ecstasy (500 ng/mL), Methadone (300 ng/mL), Methamphetamine (1000 ng/mL), Opiates (300 ng/mL & 2000 ng/mL), Oxycodone (100 ng/mL), Phencyclidine (25 ng/mL), Propoxyphene (300 ng/mL), THC (50 ng/mL) and Tricyclic Antidepressants (1000 ng/mL).

Twenty-eight different test combinations are available from single drug to 10 drugs.

3. Is this product:

a self-contained point-of-care test (requiring no equipment)?

a single-patient test for use in point-of-care equipment?

a point-of-care analyzer?

an accessory or supply used by point-of-care equipment?

other point-of-care equipment?

4. What kind of specimen/sample is tested?

Urine

5. Where is this product used? Check as many as apply:

at a community screening event

in a physician's office or outpatient clinic

at an inpatient's bedside

in an inpatient staff-access area

in a critical-care unit

elsewhere: \_\_\_\_\_

6. How long does it take to obtain test results, under ideal conditions?

5 minutes

7. How are test results made available?

A negative result displays as a color bar on the device.

8. What credentials/licensures are required to run this test?

The assay is not CLIA waived.

9. How long does it take to train the person who administers/interprets the test?

5 to 10 minutes

10. What is the typical cost:

- per test?

Per test pricing is variable (about \$4.9 per analyte) depending on the number of drugs per test. For Pricing, please contact Bio-Rad at 1-800-2BIO-RAD. For free evaluation sample, please visit [www.bio-rad.com/toxsee](http://www.bio-rad.com/toxsee).

- for initial setup (including purchase, installation, and training, as applicable)?

NONE

11. In brief, how does this product (not your company) differ most from competing products?

It is a single step assay that takes 5 minutes and uses a single use disposable dropper as opposed to other tests that require multiple pipette tip changes.