

Branan Medical Corp
Contact person: Jeff Kovalik
E-mail address: jkovalikbmc@aol.com
Phone: (949) 598-7166
Address: 10015 Muirlands Rd, Suite C
Irvine, CA 92618
Web site: www.brananmedical.com

1. What is the name of your company's point-of-care testing product or product line?
(For multiple products/lines, please complete NO MORE THAN THREE separate questionnaires.)

OratectPlusTM _____

2. What does this test detect?

Alcohol, Amphetamines, Cocaine, Marijuana, Methamphetamines/Ecstasy, Opiates, and Phencyclidine or Benzodiazepines. _____

3. Is this product:

a self-contained point-of-care test (requiring no equipment)?

a single-patient test for use in point-of-care equipment?

a point-of-care analyzer?

an accessory or supply used by point-of-care equipment?

other point-of-care equipment?

4. What kind of specimen/sample is tested?

Saliva _____

5. Where is this product used? Check as many as apply:

at a community screening event

in a physician's office or outpatient clinic

at an inpatient's bedside

in an inpatient staff-access area

in a critical-care unit

elsewhere: Hospital Lab, ER Department, Pain Management Clinic, Behavioral Health, Occupational Medicine Department.

6. How long does it take to obtain test results, under ideal conditions?

5 minutes _____

7. How are test results made available?

The test results are available on the test window of the device. The test uses a competitive immunoassay procedure involving colored colloidal gold to produce the test results. _____

8. What credentials/licensures are required to run this test?

__None_____

9. How long does it take to train the person who administers/interprets the test?

_30 minutes_____

10. What is the typical cost:

- per test? _Depends on combination of test panel and volume. _____
- for initial setup (including purchase, installation, and training, as applicable)? _No charge for training. ___

11. In brief, how does this product (not your company) differ most from competing products? _One-step oral fluid point of care test. The only oral fluid point of care test that combines drugs and alcohol on one device. _____